

PCR, Microarrays and Sequencing

Introduction to Bioinformatics Course

Niloufar Razani
Department of Computer Engineering
Sharif University of Technology



What to Read?



How to Read?

Polymerase Chain Reaction(PCR)



Microarray



Sequencing



Polymerase Chain Reaction (PCR)




Thermal Cycler

PCR Setup: What's in the Mix?

- 1) A DNA template that contains the DNA target region to amplify.
- 2) Two DNA primers that are complementary to the 3' ends of each of the sense and anti-sense strands of the DNA target.
- 3) DNA polymerase; an enzyme that polymerizes new DNA strands; heat-resistant Taq polymerase is especially common
- 4) Deoxynucleoside triphosphates, or dNTPs, the building blocks from which the DNA polymerase synthesizes a new DNA strand.

PCR Procedure



Applications


- ➤ Medical and Diagnostic
- > Forensic Applications
- > DNA sequencing can be assisted by PCR
- > Detection of Mutations
- > And so on ...

Quantitative PCR (qPCR)


- > qPCR or real-time PCR is a variation of the PCR technique.
- ➤ PCR is a qualitative technique, while qPCR is a quantitative technique.
- ➤ qPCR is based on the principle that the amount of fluorescent dye produced during a PCR reaction is proportional to the amount of DNA that is amplified.

SYBR Green is used in qPCR in the following way:


- 1. The SYBR Green dye is added to the PCR mixture.
- 2. The PCR reaction is then performed as usual.
- 3. The fluorescence of the SYBR Green dye is measured at each cycle of the PCR reaction.
- 4. The amount of DNA that has been amplified is determined by measuring the increase in fluorescence over the course of the PCR reaction.



qPCR Curve




qPCR Curve




RT-qPCR

- Reverse transcription qPCR (RT-qPCR) is a technique that combines reverse transcription (RT) and quantitative PCR (qPCR) to detect and quantify RNA.
- ➤ RT is a process by which RNA is converted into complementary DNA (cDNA).
- Common applications of RT-qPCR include gene expression analysis, mRNA quantification, RNA virus detection and so on ..




Microarray

- ➤ A DNA Microarray is a collection of microscopic DNA spots attached to a solid surface.
- > DNA microarrays are used to measure the expression levels of large numbers of genes simultaneously.
- Each DNA spot contains a specific DNA sequence, known as probes (or reporters or oligos).
- These can be a short section of a gene or other DNA element that are used to hybridize a cDNA.





Microarray



Sequencing - Sanger


- Sanger sequencing is based on the use of dideoxynucleotides (ddNTPs).
- ➤ ddNTPs are similar to normal nucleotides(dNTPs), but they are missing a hydroxyl group on the 3' carbon atom.
- This prevents the DNA polymerase enzyme from adding another nucleotide to the chain, and it therefore terminates the chain.





Sanger Sequencing Steps

- 1. Prepare the DNA template. The DNA template is the DNA sequence that you want to sequence. It can be a single gene, a whole genome, or any other region of DNA.
- 2. Prime the DNA template. A primer is a short piece of DNA that binds to the beginning of the DNA template. The primer is used to start the DNA polymerase reaction.
- 3. Add the DNA polymerase enzyme, nucleotides, and ddNTPs. The DNA polymerase enzyme is responsible for adding nucleotides to the DNA chain. The nucleotides are the building blocks of DNA, and the ddNTPs are the chain-terminating nucleotides.
- 4. Perform the DNA polymerization reaction. The DNA polymerase enzyme will add nucleotides to the DNA chain until it encounters a ddNTP. At that point, the chain will terminate.
- 5. Separate the DNA fragments by size. The DNA fragments of different sizes are separated using a technique called gel electrophoresis.
- 6. Detect the fluorescent labels. The ddNTPs are labeled with different fluorescent dyes. This allows the DNA fragments to be detected and visualized under ultraviolet light.



Sequencing - Sanger



Next Generation Sequencing

First generation

Second generation (next generation sequencing)

Third generation

















Sanger sequencing Maxam and Gilbert Sanger chain termination

Infer nucleotide identity using dNTPs, then visualize with electrophoresis

500-1,000 bp fragments

454, Solexa, Ion Torrent, Illumina

High throughput from the parallelization of sequencing reactions

~50-500 bp fragments

PacBio Oxford Nanopore

Sequence native DNA in real time with single-molecule resolution

Tens of kb fragments, on average

Short-read sequencing

Long-read sequencing

NGS Cost

